

Objective

To characterize the microbial killing efficacy of sterilants utilizing ATS.

Methods

To characterize the microbial killing efficacy of the sterilants to be evaluated under worst-case conditions, 10 uL of ATS containing $\sim 10^8$ cfu/ml of the test microorganism is inoculated inside a surrogate lumen carrier (3mm internal diameter, 2 cm long), and dried at room temperature overnight. The inoculated test segment is then assembled into a lumen of 125cm total length and then processed through the sterilization process. The test carriers consist of PVC plastic tubing cut into 2 cm lengths or PTFE endoscopy biopsy channel tubing cut into 2 cm lengths.

Sterility Testing

The central lumen segment is aseptically removed and placed in a tube containing 2 mLs of sterile tryptic soy broth containing 10% fetal bovine serum. This segment is incubated at 35°C for 5 days for *Enterococcus faecalis*; 55°C for 5 days for *Geobacillus stearothermophilus*; and 30°C for 10 days for *M. chelonae*. Because turbidity may be hard to detect for *M. chelonae*, blind subcultures are performed after 10 days of incubation by sub-culturing the broth to blood agar (BA) and incubating the plates at 30°C for 5 days.

To quantify the level of residual viable test organisms post-treatment, the entral segment of the duplicate set of lumen test carriers is aseptically disassembled and placed into a sterile tube containing 2 mLs of sterile tryptic soy broth containing 10% fetal bovine serum. The tube is sonicated in a Branson 1220 sonicator bath for 2 x 5 secs, and then mixed on an S/P Multitube vortexer (American Dade, Miami, FL) on setting "2" for 10 minutes. Serial 1:10 dilutions are prepared in sterile tryptic soy broth and 0.1 mLs plated on tryptic soy agar plates (BA plates for *M. chelonae*) using the spread plate technique. The inoculated plates are incubated at 35°C for 24 - 48 hours for *E. faecalis*, 55°C for 24 - 48 hours for *G. stearothermophilus* and 30°C for 48 - 72 hours for *M. chelonae*. Colonies are enumerated and the viable count per test carrier is calculated.