



SAFETY DATA SHEET

SDS FOR BactiQuik™ MICROBIAL SUSPENSIONS

(Biosafety Level 1)

BactiQuik™ Microbial Suspensions are not hazardous. However, as living microorganisms they are potentially hazardous.

Description

BactiQuik™ Microbial Suspensions consist of microorganisms suspended in liquid medium. Suspensions are shipped in liquid form at refrigerated temperatures.

SECTION 1

Hazardous Ingredients

Suspensions do not contain any hazardous ingredients.

SECTION 2

Physical Data

Liquid suspensions shipped in glass vials.

SECTION 3

Health Hazards

This microbial suspension is not known to cause disease in healthy human adults or animals.

SECTION 4

Fire and Explosion

Not applicable.

SECTION 5

Reactivity Data

Stable. Hazardous reactivity will not occur.

SECTION 6

Method of Disposal

Spill: Contain the spill and decontaminate spill area using suitable disinfectants such as chlorine bleach or 70% ethyl or isopropyl alcohol. Product Disposal: Dispose any unused portion of the product by autoclaving at 121° C for a minimum of 20 minutes.

Follow all Federal, State, and local regulations.

SECTION 7

Special Protection Information

For Biosafety Level 1 Microbial Suspensions

Handle as potentially biohazardous material under at least Biosafety Level 1 containment as defined by the Center for Disease Control Office of Health and Safety.

SECTION 8

Precautions and Comments

Healthmark Ind. recommends that qualified microbiologists using appropriate safety procedures handle BactiQuik™ Microbial Suspensions. Further information regarding laboratory safety procedures can be obtained from the Center for Disease Control Office of Health Safety's web site at www.cdc.gov/od/ohs/biosfty/bmb14toc.htm.

THE ABOVE INFORMATION IS CORRECT TO THE BEST OF OUR KNOWLEDGE. ALL SUSPENSIONS MAY PRESENT UNKNOWN HAZARDS AND SHOULD BE USED WITH CAUTION. THE USER SHOULD MAKE INDEPENDENT DECISIONS REGARDING THE COMPLETENESS OF THE INFORMATION BASED ON ALL SOURCES AVAILABLE. HEALTHMARK SHALL NOT BE HELD LIABLE FOR ANY DAMAGE RESULTING FROM HANDLING OR CONTACT WITH THE ABOVE PRODUCT.